




INTERFINE CHEMICALS

ANALYTICAL SCIENCES

BIOCHROMATOGRAPHY

LIFE SCIENCES

Tél. +33 (0)4 70 03 73 01 ■ interfine@interchim.fr

Technical Data Sheet

PURIFLASH® 4250₂₅₀

250 ml/min @ up to 250 bar


Premium preparative LC

Lab Scale

- ✓ Flash columns up to F800 format
- ✓ prep-LC columns up to 80mm id - 5µm

puriFlash® 4250 is a Premium, easy to use prep LC system. Perfect for routine & sophisticated purifications with up to 5µm & up to 80mm id prep LC columns.

Interchim® Software ver. 5.0



One customizable screen to manage the whole process to create the method, to launch the purification & to control the process during the purification up to the report. User Profile definition: Customization of the GUI and its functions User profile management and password Definition of the reporting format, etc... Unique features: TLC to Flash transposition, "Best chromatogram", Automatic gradient optimization, Intelligent Continuous Mode, Dual triggered collection, Timed auto shut-off, Air purge...etc...

P/N: KV6090



Options...

Detector

- UV 190 - 840 nm - (GO9630)
- puriFlash® ELSD - (F10320)
- puriFlash® Mass Spec - (JO1020)
- RI - (FSQ600)

Valves

- Column selection - (LV2650)

Column Holder

- Stand Alone - (PF4530)

Fraction collector

- Drainage system of the tray collector - (DZ7790)
- Fume enclosure - (PF4350)
- Funnel rack set - (FJ4330)

Unit Control

- Bar code reader - (GO2380)
- CCD camera - (GO2390)

Safety solvent

- Solvent level sensor Global - (FSP720) or Solvent level sensor Separated - (FSP730)
- Safety solvent bottle caps - (DV2760)
- Safety waste cap - (IO6930 or JO4500)

Open Access Lab

- 10 columns sequential - (JO3750)

INTERFINE CHEMICALS

ITM-2012/08/09



interchim®

® 211 bis, avenue JF Kennedy
BP 1140 - 03100 Montluçon
Fax 33 (0)4 70 03 82 60

■ Hotline +33 (0)4 70 03 73 01 • interfine@interchim.fr ■



