



CMA 120 Bowl with Food and Water Containers, CMA 120 Sviwel Assembly and Balance Arm.

- External containers do not disturb the animal's rotation
- Easy to disconnect for cleaning

# CMA 120 Bowl with Food & Water Containers

for use with a CMA 120 Freely Moving System

The CMA 120 Bowl with Food and Water Containers is used in studies where a freely moving animal will be confined for longer periods of time. The additions of these containers allow the animal to feed and drink ad libitum.

The food container and water bottle are arranged on the outside of the bowl so as not to disturb the movement of the tethered animal when inside the bowl. Both containers are easily removed for cleaning and refilling.

## Technical information

### Bowl

Height: 360 mm  
Diameter: 400 mm  
Material: Perspex

### Water Bottle

Volume 150 mL

### Food Container

Diameter ID: 5 cm, OD: 6 cm

## Ordering information

CMA 120 Bowl with Food and Water Containers Ref. No. 8309672

