

## Basket Shafts



QLA Basket Shafts are manufactured from electropolished 316 Stainless Steel. Clip on style and O-ring style shafts are available for most dissolution instruments. O-ring style shafts comply with Original Equipment Manufacturer's dimensional specifications and are guaranteed to properly interface with the designated dissolution basket. Clip on style shafts utilize a proprietary, spring tempered clip, which retains the original gripping tension far longer than conventional shafts. Clip on style shafts comply with USP (711) and Original Equipment Manufacturer's dimensional specifications. Both clip on style and O-ring style shafts are precision manufactured to limit rotational wobble to levels below USP (711) specifications.