



## Chromeleon 7



### [Latest News about Chromeleon 7](#)

The Dionex Chromeleon 7 [Chromatography Data System](#) delivers everything you need to get the most from UHPLC. It lets you create your analytical experiments faster, process all your data in the shortest possible time, and instantly create your results. All this is delivered through rich, intelligent functionality in combination with Operational Simplicity. Chromeleon 7 provides solutions for both workstation and enterprise customers. It can control a single LC, IC, or GC and be easily extended to control hundreds of these instruments. Its unique XVault(TM) architecture ensures that data is always secure and instruments can always be used – even in the event of network outages. Chromeleon 7 gets you from samples to results quickly and easily. Combined with UltiMate 3000 systems this will boost your overall lab productivity.

**Chromeleon 7 is Simply Intelligent.**

## Pełna specyfikacja

### Streamline your workflow

The Chromeleon 7 eWorkflow framework minimizes operator tasks between sample receipt and generation of the final results, dramatically increasing lab productivity. Choose an eWorkflow, select an available instrument, enter the number of samples and start the analysis. The eWorkflow automatically creates the sequence, starts the run, and ensures that data are processed and reported correctly, minimizing user error and reducing review time.

### Gain Instrument Control

All your UltiMate 3000 systems and other systems can be controlled through customizable ePanels, providing access to all instrument control features; Entering basic instrument parameters, like flow, temperature and detector wavelengths or more advanced features, like monitoring the instruments wellness parameters. With all your instruments in the lab always readily available in just one click instrument control has never been easier.

### Speed Up Data Analysis

Peak detection can be a major decision area in data analysis. For the challenge of unresolved peaks, the unique and intuitive SmartPeaks Integration Assistant gives you the baseline you want. Simply select a region of the chromatogram and choose the integration treatment you want. No training is required, so everyone can immediately integrate complex chromatograms. SmartPeaks saves time and effort and moves you quickly to your integration goals.

### Chromeleon Tips and Tricks

Dionex offers short clips that provide tips on how to get the most from Chromeleon. Some great examples are

Creating a Summary Trend Report

How to find a Chromeleon Formula

Changing the Channel in a Miniplot



Converting an HPLC method to UHPLC