

CMA 110 Liquid Switch

- instantly switch between syringes
- prevents introduction of air bubble
- no interruption of flow

The CMA 110 Liquid Switch permits manual switching between up to three perfusion lines (syringes) and a microdialysis probe. This makes it possible to change different solutions instantaneously without any risk of introducing air bubbles into the microdialysis probe.



Ordering information

Ordering information for CMA 110

Ref. No.

830 8200 CMA/110 Liquid Switch. Tubing Kit included

Accessories:

340 8550 CMA 110 and 111 Tubing Kit