

## CMA 110 Liquid Switch

- instantly switch between syringes
- prevents introduction of air bubble
- no interruption of flow

The CMA 110 Liquid Switch permits manual switching between up to three perfusion lines (syringes) and a microdialysis probe. This makes it possible to change different solutions instantaneously without any risk of introducing air bubbles into the microdialysis probe.



---

### Ordering information

---

#### Ordering information for CMA 110

---

##### Ref. No.

830 8200 CMA/110 Liquid Switch. Tubing Kit included

##### Accessories:

340 8550 CMA 110 and 111 Tubing Kit