

## CMA 111 Syringe Selector



- instantly switch between syringes
- prevent introduction of air bubble
- no interruption of flow

The CMA 111 Syringe Selector allows automated or manual switching between one of three perfusion lines (syringes) and a microdialysis probe. This makes it possible to change different solutions instantaneously without any risk of introducing air bubbles into the microdialysis probe.

The CMA 111 Syringe Selector is motor-driven and powered from the CMA 400 Microinjection Pump.

---

### Ordering information

---

### Ordering information for CMA 111

---

#### Ref. No.

830 8210 CMA 111 Syringe Selector. Tubing Kit included

#### Accessories:

340 8550 CMA 110 and 111 Tubing Kit  
800 2030 CMA 400 Adapter for CMA 111