

Membrane nitrogen generator – MNG series



The Membrane nitrogen generator has been developed to meet specific requirements in terms of flow, purity and pressure for LC-MS applications. It can also be used for the evaporation of solvents in samples being analysed. The simple high efficiency membrane technology allows the separation of nitrogen from the other components of the compressed air inlet.

Specyfikacja

RANGE	REFERENCE	OUTLET FLOW RATE	PURITY	OUTLET PRESSURE
MNG Series: high efficiency membrane technology	MNG8	8 l/min	> 99%	7bar
	MNG40	40 l/min	> 99%	7bar
	MNG80	80 l/min	> 99%	7bar
	MNG120	120 l/min	> 99%	7bar