

## Microsyringes



Microsyringes 1, 2.5, 5 and 10 mL are available. The Microsyringes are glass syringes with piston stroke 60 mm. Syringe Needles are pretreated to reduce degradation of labile molecules, such as catecholamines and their metabolites during a microdialysis experiment.



### Ordering information

### Ordering information for Microsyringes

Ref. No.

830 9020 Microsyringe 1 mL glass

830 9021 Microsyringe 2.5 mL glass

830 9022 Microsyringe 5 mL glass

830 9023 Microsyringe 10 mL glass \*

743 1083 Syringe Needle, 5/pkg

\* for use with CMA 400 Syringe Pump