



## Nitrogen / Zero Air



The combine N<sub>2</sub>/Air generator eliminates the need for costly, inconvenient high pressure cylinders in the laboratory. Including an integral oil free air compressor in option, the combine nitrogen/air generators deliver a continuous stream of pure nitrogen gas > 99,9995% and grade purified air with HCs content <0,1ppm, ideal for applications, such as: GC carrier gas, make up and GC-FID.

### Specifications

SPECIFICATIONS	ANG1 ZANG1*	ANG3 ZANG3*
Max. N <sub>2</sub> (> 99,9995%) flow rate	1L/min	3L/min
HCs < 0.05 ppm for N <sub>2</sub> gas	- yes	- yes
Purified Air (HCs < 0,01 ppm, H <sub>2</sub> O dp < - 55°C) flow rate	1.5L/min	3L/min
Outlet pressure	5,5 bar	