

## Perfusion Fluid



### PERFUSION FLUID CNS CONTENTS

Perfusion Fluid is a sterile, isotonic fluid especially developed for Microdialysis use.

It is conveniently packed in glass ampoules, 5 mL, ready to use.

Perfusion Fluid CNS for studies in brain tissue.

### PERFUSION FLUID T1 CONTENTS

Perfusion Fluid is a sterile, isotonic fluid especially developed for Microdialysis use.



---

## Specyfikacja

### PERFUSION FLUID CNS CONTENTS

NaCl	147mmol/L
KC	2.7 mmol/L
CaCl 2	1.2 mmol/L
MgCl 2	0.85 mmol/L

### PERFUSION FLUID T1 CONTENTS

NaCl	147mmol/L
KC	4 mmol/L
CaCl 2	2.3 mmol/L
MgCl 2	0.85 mmol/L