

Podstawowy HPLC



Thermo SCIENTIFIC

The UltiMate 3000 [Basic Automated LC](#) system is optimized for reliability and ease-of-use with routine LC applications, and also offers full UHPLC compatibility. The ACC-3000 Autosampler Column Compartment is at the heart of the system. Its unique instrument design combines a rugged sample injection principle with a powerful column oven.

- Future proof yet cost-effective alternative to standard systems
- 620 bar maximum pressure and 100 Hz data rate for UHPLC compatibility
- Minimum number of wear parts and critical connections
- Industry leading range of detectors
- Modular versatility at compact system size