



## PSA nitrogen generators – Maximus Lite series



The nitrogen generators serie Maximus Lite offer to you a unique, innovative solution to nitrogen gas supply. A reliable, secure source of nitrogen can be produced from your existing compressed air supply, eliminating the need for liquid nitrogen or high pressure gas cylinders. Conceived for a continuous operation, the MAXIMUS LITE nitrogen generator can produce flow rates from a few liters per minutes to over NL/min at purities in oxygen content residual from 10 ppm to 3%.

### Specyfikacja

RANGE	REFERENCE	OUTLET FLOW RATE	PURITY	OUTLET PRESSURE
Maximus Lite series	MNG102L	10 L/min	> 99,999%	5-14bar
PSA technology	MNG104L	20 L/min	> 99,999%	5-14bar