

SPE columns



Thermo
SCIENTIFIC

Solid Phase Extraction Columns

A range of versatile, high quality SPE solutions for rapid, effective and targeted sample preparation.

The highly reproducible and efficient phases provide consistently high recoveries, free from contaminants and impurities.

SOLA | HyperSEPSERVO+ | HyperSEP

Version

SOLA

Thermo Scientific™ SOLA introduces a revolutionary change in SPE technology providing unparalleled performance that results in greater confidence in analytical results, first time success and ultimately cost and time saving. In today's high throughput bioanalytical and clinical research laboratory environment where reproducibility, certainty of results and cost saving is a fundamental requirement, SOLA products are an indispensable tool to provide confidence and first time/every time success in the analytical process.

Columns: HRP SCX, SAX.

HyperSEPSERVO+ and SERVO

Both ranges offer high levels of reproducibility and cleanliness of extract for dealing with difficult sample matrices where extra confidence in your results is required. Servo+ provides the additional benefit of greater selectivity, higher loading capacities and robustness.

HyperSEP Servo+ and Servo offer dedicated solutions for the analysis of amphetamines, cocaine, opiates, PCP, and THC as well as total screening.

Columns: Total, Total A, Total B, THC, Opiates, AMP, PCP, Cocaine.

[Brochure](#)

HyperSEP

Over the last few years, the requirement for higher sensitivity, selectivity, accuracy, precision and sample throughput has increased significantly. This is due to reduced sample volumes, greater drug efficacy, and greater awareness of toxicity levels of pollutants.

Columns: Retain PEP, Retain-CX, Retain-AX, Hypercarb, C18, C8, Phenyl, Silica, Florisil, Cyano, Diol, Aminopropyl, SAX, SCX,



Verify-CX, Verify-AX.

[Brochure](#)