

Stereotaxic Instruments for Mouse



Just for Mouse™ Stereotaxic Instruments

THE JUST FOR MOUSE™ STEREOTAXIC INSTRUMENTS is ideal for researchers working with mice and other small rodents. Precision alignment when using the Just for Mouse™ Stereotaxic ensures accurate placement of electrodes, micropipettes, cannula and other devices.

Advanced digital and motorized versions include a zeroing function allowing the user to target specific coordinates in the brain for injection, implantation, etc. Calculation of distance measurements is thus no longer necessary.

The Just for Mouse™ Stereotaxic Equipment is available in manual, digital and motorized versions.

Ordering information

Ordering information for Just for Mouse™ Stereotaxic Instruments

Ref. No.

CMA 72-9564 Just for Mouse™ Stereotaxic Model 51730

CMA 72-9565 Dual Just for Mouse™ Stereotaxic Model 51733

CMA 72-9563 Digital Just for Mouse™ Stereotaxic Model 51730D

CMA 72-9566 Digital Dual Just for Mouse™ Stereotaxic Model 51733D

CMA 72-9567 Motorized Just for Mouse™ Stereotaxic Model 51730M