

## Szafki na kolumny



### Szafki na kolumny chromatograficzne - Columnstore

The ColStor storage system is simple, compact and stops columns getting lost or damaged. Poor column storage can lead to reduced column performance and decreased lifetime.

- Stores 30, 60, 90 or 120 HPLC columns - in secure, moulded drawer inserts
  - Optional Fittings Tray - Moulded 6 position tray for easy storage of HPLC tools and fittings
  - Stackable — Unique stacking kit allows you to increase the number of columns you can store at any time.
  - Holds columns up to 30cm long - will also take 5cm, 12.5cm, 15cm and 25cm!
  - 100% Steel Cabinet — Quality, solid steel construction with chrome-plated D-Ring handles
  - Eliminates Bulky Column Boxes — Columns held securely in moulded drawer inserts
  - Holds all end-fitting types — Narrow and wide end-fittings held in flexible foam
  - Versatile — Store up to 12 x 12.5cm columns in a single drawer, or a combination of short columns, guard columns and analytical columns
- The Labhut Column Store is a cost-effective, expandable solution to properly indexing and storing your HPLC columns

