



TOC AIR GENERATOR – Serie TOC



F-DGS's TOC Air generators are some of the smallest and most efficient Hydrocarbon-free air and CO₂ purifier systems for laboratory use. Utilising a heated catalyst technology which is combine with a CO₂ removal air dryer to remove THC including methane to less than 0.05 ppm and CO₂ less than 1 ppm. This TOC air station is controlled using the latest in HMI touch screen technology to display the process in real time.

Specyfikacja

RANGE	REFERENCE FLOW RATE		PURITY	OUTLET PRESSURE
TOC series: heated catalyst technology	TOC -1	1500 ml/min	HC < 0,05 ppm	up to 6.0 bar
	TOC-3	3000 ml/min	CO < 0,05 ppm CO ₂ < 1 ppm H ₂ O dp< -70°C	