

ULTRA ZERO AIR GENERATOR - Serie ZAG



F-DGS's Ultra Zero Air Generators are some of the most efficient Hydrocarbon-free air purifier systems for laboratory use. Utilising the PSA technology to remove water vapour and the heated catalyst technology to remove THC including methane to less than 0.05 ppm. This Ultra Zero Air station DO NOT requires external compressed air supply. The Generator is controlled using the latest in HMI touch screen technology to display the process in real time.

Specyfikacja

| RANGE | REFERENCE FLOW RATE | PURITY | OUTLET PRESSURE |
|-------------------------------|---------------------|---|-----------------|
| ZAG series: PSA technology | ZAG-1 | 1500 ml/min | up to 6.0 bar |
| | ZAG-3 | 3000 ml/min | |
| | ZAG-6 | 6000 ml/min | |
| | ZAG-10 | 10 L/min | |
| | ZAG-20 | 20 L/min | |
| | ZAG-25 | 25 L/min | |
| | ZAG-35 | 35 L/min | |
| | | CH ₄ , CO, NO _x < 0,05 ppm | |
| | | CO ₂ < 5 ppm | |
| | | H ₂ O dp < -50°C | |