

ULTRA ZERO AIR GENERATOR – Serie ZAG



F-DGS's Ultra Zero Air Generators are some of the most efficient Hydrocarbon-free air purifier systems for laboratory use. Utilising the PSA technology to remove water vapour and the heated catalyst technology to remove THC including methane to less than 0.05 ppm. This Ultra Zero Air station DO NOT requires external compressed air supply. The Generator is controlled using the latest in HMI touch screen technology to display the process in real time.

Specyfikacja

RANGE	REFERENCE FLOW RATE	PURITY	OUTLET PRESSURE
ZAG series: PSA technology	ZAG-1	1500 ml/min	CH ₄ , CO, NO _x < 0,05 ppm CO ₂ < 5 ppm H ₂ O dp< -50°C up to 6.0 bar
	ZAG-3	3000 ml/min	
	ZAG-6	6000 ml/min	
	ZAG-10	10 L/min	
	ZAG-20	20 L/min	
	ZAG-25	25 L/min	
	ZAG-35	35 L/min	