

ZERO AIR GENERATOR – Serie ZA



The Zero Air generators produce dry and hydrocarbon-free air, using air from oil-free compressors, thus avoiding the need to use traditional bottles that are often complex to change. Designed with safety and convenience in mind, this system will eliminate the need for inconvenient high pressure gas cylinders. This complete turnkey system is engineered with the highest quality components, is easy to install and requires only minimal annual maintenance.

Specyfikacja

RANGE	REFERENCE	FLOW RATE	PURITY	OUTLET PRESSURE
ZA series:	ZA-1500	1500 ml/min	HCs < 0,05 ppm CO < 0,05 ppm	up to 6.5 bar
	ZA-3000	3000 ml/min		
	ZA-6000	6000 ml/min		
	ZA-15000	15 L/min		
	ZA-30000	30 L/min		